



PRODUCT CATALOGUE



Model Production

TOPDEB Debubbler	5
Top Vacudeb	
Top Sep	
Top Blockout Wax	
Top shield&Shine Spray	6
Top Smart Sep	
Top Die Hardener	
Top Die Hardener - Pen Applicator	7
Top Die Sealer	

Die Preparation

Top PC Spacer	
Top PC Spacer - PCS ISO	9-10
Top PC Spacer - KITS	
Top Die Spacer	
Top Die Lube	11
Top Die Lube - Pen Applicator	
Top Rubbsep Top Blockout Wax	12

Ceramic Products

Top Ceraliq Top Cera Glazer Top Marker	15
Top Cerafix Gel Top Cera Contour	16
Top Die Hardener Top Die Hardener - Pen Applicator Top Die Sealer	17
Top Cera Sep Top Cera Sep - Pen Applicator Top Cera Sep Spray	18
Top PC Spacer Top PC Modeling Wax	19
Ceramic Product Kits	20

DST-FLEXI DENTURE MATERIALS

DSTVITOflexi	23
DSTVITOacry	24
DSTVITOcetal	25
DSTVITOkeep	26
DSTVITOpplast	27-28
DSTVITOlux-f	29
DST CLEAN & RENEW DST LUSTERLIQ DST LUBRILIQ	30
DST ISOLIQ F DST CEM F DST FUSE FLEX	31

Denture

Top Natur Glaze VITO CRYL Seal -IQ®	34
Top Air Opaque Top Air Opaque - Top Paqubond Top Air Opaque - KITS	35 - 36
Top Airstain Glaze Top Airstain Glaze - KITS	37 - 38
Top Pumiliq Top Gyps Rid Top Blockout Wax Top Bite Rims	39
Top Sticky Wax Top Spec Sticky Wax Top Isoliq A	40

CAD/CAM

Top Scan Wax Bari x-ray	43
----------------------------	----

WAXES

Dipping Waxes - Top Dip Wax S Top Dip Wax M Top Dip Wax H	45
Crown and Bridge Waxes - Top Modeling Wax Top PC Modeling Wax Top Sculpturing Wax	46
Special Waxes - Top Diagnostic Wax Top Cervical Wax Top Milling Wax Top Scan wax Top Blockout Wax	47- 48
Bite Registration Waxes - Top Bite Rims	48
Sticky Wax - Top Sticky Wax Top Spec Sticky Wax	49

CASTING/SOLDERING

Top Wax Finish Uni-Wax Adhesive Top Vacudeb	51
---	----

DS Technology LTD

INNOVATION THAT MAKES THE DIFFERENCE

Developing a range of dental solutions requires talent, experience, and an eye for detail. Headed by a skilled dental technician with 22 years of hands-on experience, DS Technology makes the very best products available to dental technicians and dentists.

At **DS Technology**, we have carefully manufactured a comprehensive line of competitively priced, quality dental materials designed especially for dental technicians and dentists. Each material improves work efficiency, enhances final results and increases endurance, resulting in a highly reliable application.

We've been there on the front lines, in busy practices with anxious patients and overloaded professionals. We know what it's like to discover a process that shaves off precious minutes, stabilizes faster and provides reliable, high quality results. Our specialized products and technologies have made procedures more efficient, accurate and practical for dental professionals, as well as more comfortable and with better fit and finish for their patients.





Model Production

Top Deb - Debubblizer | Top Vacudeb | Top Sep
Top Die Hardener | Top Die Sealer | Top Blockout Wax
Top Smart Sep | Top shield & Shine Spray



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

TOPDEB Debubbler



When pouring into molds for casting dental models, surface tension can push back the casting material and form imperfections and air bubbles. Our Top-Deb Debubbler prevents this by stabilizing the surface of the mold, allowing material to flow into the thinnest areas and adhere completely to the mold, allowing accurate, highly detailed and dense castings. Can also be used in treating stored impressions and wax patterns, to aid in removal of dust particles and oily residues which accumulate over time.

- | Universal application, including silicone, rubber base or alginate impression material
- | Transparent, blue tint.
- | Includes refillable misting spray container for easy application.

Part num	Color	Amount
26-400/250	Blue - Transparent	8 oz
26-400/500	Blue - Transparent	16 oz
26-400/1000	Blue - Transparent	32 oz

Top Vacudeb



A surface tension reducer, Top Vacudeb allows an attachment to achieve close adaption and maximum contact with all surfaces of a wax pattern and with all impressions prior to investing, without trapping air, forming air bubbles or voids on the surface of models and castings.

By managing surface tension, an improved flow rate of investment material and strict adherence to the impression ensures that the casting is reproduced perfectly and accurately.

- | May be used under vacuum.
- | Fast drying increases productivity.
- | Eliminates air entrapments, bubbles and voids.
- | Facilitates fluid flow to produce sharper detail.
- | Supports silicone, wax, rubber and alginate patterns.
- | Transparent pink tone.

Part num	Color	Amount
26-220/250	Pink- Transparent	8 oz
26-220/500	Pink- Transparent	16 oz
26-220/1000	Pink- Transparent	32 oz

Top Sep



A non-oily anti-sticking separating agent, and dries on contact. Since it is alcohol-based, Top Sep quickly dries to a hardened finish on the treated surface, allowing immediate pouring of investment into dies or molds. There is no film leave-behind.

- | Environmentally safe; contains no silicones or chlorofluorocarbons.
- | Soaks into stone or plaster to seal pores and assure smooth finish of investment.
- | Hardens and smooths surface of stone or plaster, and makes it water-proof.
- | Quick cleanup if sprayed on trays, vessels, flasks, sinks, etc.
- | Leaves no film.

Part num	Color	Amount
26-200-250	Yellow Transparent	8 oz
26-200-500	Yellow Transparent	6 oz
26-200-1000	Yellow Transparent	32 oz

Top Smart Sep

TOP SMART SEP specially formulated for separating casting and modeling waxes from porcelain, stone, composite, cement, acrylic resins and metal.



- | Universal application.
- | Easy to apply.
- | Restorations exhibit smoother and cleaner surfaces.

Part num	Color	Amount
26-201/8	Transparent	8oz
26-201/16	Transparent	16oz
26-201/32	Transparent	32oz

Top Shield & Shine

High gloss coating to all gypsum models - from orthodontic to the practical models



- | Fast drying
- | Easy to apply
- | Get the professional shine luster and hardens
- | Prevents staining & chipping of models
- | Provides a smooth, glossy stone models
- | Ideal from orthodontic models
- | Professional quality look when presenting models - improves case presentations.
- | Works on all dental stones & plasters
- | Thickness of less than one micron
- | Accuracy & quality are maintained
- | Not create occlusal inter-ference.

Part num	Color	Amount
26-925/4	Transparent	4oz
26-925/8	Transparent	8oz
26-925/16	Transparent	16oz
26-925/32	Transparent	32oz

Top Blockout Wax

A very hard, high temperature wax, our Top Blockout Wax is designed to fill undercuts and process induced cavities, cover sharp edges and smooth out rough surfaces of dies prior to molding attachments. Excellent carvability and very good scraping properties.



- | Fast and secure application to undercuts on dies for molding of comfortable attachments.
- | Adapts to die contours while soothing.
- | Does not burn without residue.
- | Will not coalesce with dipping wax, cervical wax or other modeling waxes due to very high melting point.
- | Does not discolor plaster model after boil-out.
- | For blocking out bite guard splints and dies.
- | 105 °C melting point.

Part num	Color	Amount
26-500	White Transparent	42/45g
26-501	Pink Transparent	42/45g



Top Die Hardener

Our Die Hardener is an absorbing resin that seals, hardens and strengthens stone models and dies. As it is absorbed into the stone surface, it permeates into all of the fissures and microscopic pores, bonding the entire stone by fusing the particles together. It resists chipping by leaving a zero micron infusion to the surface layer bonded to the resin permeating the stone, assuring an accurate and durable mold.



- | Separating dies for easy removal from investments
- | Lasts longer than other sprays
- | Pleasant smell
- | Environment friendly
- | Quick and easy wax clean up
- | Allows immediate pouring of investment
- | Spray form
- | Contains no silicons or chlorofluorocarbons

Part num	Color	Amount
26-930	Transparent	1 oz
Thinner		
26-911	Transparent	1 oz

Special thinner for setting preferred consistency Especially adapted to the texture of the TOP DIE HARDENER & SEALER

Top Die Hardener Pen Applicator

- | Reliable separating efficiency.
- | Comfortable and Easy application due to pen shape.
- | Economical and on-the-spot application;
- | Fast drying in only seconds. Hardens, seals, smooths
- | Fills open micro pores on the plaster surface
- | Abrasion-resistant surface. Dirt&water-repellent
- | Resists chipping and cracking.

Part num	Color	Amount
26-931	Transparent	1/2 oz
Refill		
26-932	Transparent	100ml



Top Die Sealer

Our Top Die Sealer is brush applied to easily seal and water-proof stone die and model surfaces. It protects margins and reduces abrasion during operations, protecting the die from chipping, breaking or sticking. It also can provide an excellent base layer for our Top Cera Sep separating substance between porcelain and cast. It separates porcelain from model material on gingival surface of pontics, in both marginal porcelain regions as well as adjacent contact zones.



- | Protects preparation borders of crowns and nlays.
- | Prepares porcelain shoulder layers.
- | Adds negligible film thickness.
- | Thin consistency allows penetration into stone for firm seating.
- | Transparent surface becomes hard and water-proof.
- | Requires only one application, although additional thin layers can be applied after first coat is dry.

Part num	Color	Amount
26-920	Transparent	1 oz
Thinner		
26-911	Transparent	1 oz





Die Preparation

Top PC Spacer | PCS ISO | Top Die Spacer | Top Die Lube
Top Rubbsep | Top Blockout Wax



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top PC Spacer

FOR PRESSABLE CERAMIC TECHNIQUE

For the creation of highly aesthetic all-ceramic restorations, including inlays and onlays, veneers and partial crowns. Enhances the visual perfection by emulating the natural shade. Its natural tint provides optimum compatibility with an accurate shade matching the patient's dentition to perfect results

Applies quickly and evenly , Dries quickly , Smooth surface, Simulating accurate die cement spacing.

ATTENTION: Use only with PCS ISO - special separating agent

The use of incompatible isolating agents could provoke unforeseen reactions

Available in : 9 IPS® Natural Die Material® shades * and milk shade, 1/2oz

*IPS Natural Die Material® is a registered trademark of Ivoclar Vivadent.



Part no	52-2801	52-2802	52-2803	52-2804	52-2805	52-2806	52-2807	52-2808	52-2809	52-2810
Color	PCS 1	PCS 2	PCS 3	PCS 4	PCS 5	PCS 6	PCS 7	PCS 8	PCS 9	PCS Milk
Amount	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz

Top PCS Thinner - Part no -52-285, Color - Transparent, Amount - 1/2 oz

PCS ISO

PCS ISO is specially formulated as a compatible preparation with Top PC Spacer and recommended to be used together. The use of incompatible isolating agents could provoke unforeseen reactions. Forms a very thin layer which prevents surface falsification. PCS ISO is formulated from environmentally safe materials, contains no solvents and will not attack die varnishes or waxes. It won't evaporate either.

ADVANTAGE IS a optimum compatibility and NOT damage the - TOP PC SPACER- TOOTH COLOURED DIE SPACERS



- | Facilitates the removal of wax patterns .
- | Prevents the die material from sticking to the ceramic material.
- | Ensures precise and comfortable fit of wax couplings.
- | Developed for use with our Top PC Spacer for speedier and simple application.
- | Low surface tension

Part num	Color	Amount
52-284	transparent	1/2 oz

Top PC Spacer -KITS



COMPARISON CHART

Top PC Spacer Shades	PCS1	PCS2	PCS3	PCS4	PCS5	PCS6	PCS7	PCS8	PCS9	PCS MILK
IPS® Die Natural Material Shades	ND1	ND2	ND3	ND4	ND5	ND6	ND7	ND8	ND9	
Vita®	A1,A2	A1,A2 A3, B2	B3,B4	A3.5	A4		C1,C2, C3, D1, D2, D3			

COMPLETE KIT

Part num - 52-280
 All 10 Natural Shades in 1/2 oz bottle: PCS 1, PCS 2, PCS 3, PCS 4, PCS 5, PCS 6, PCS 7, PCS8, PCS 9, PCS MILK
 PCS Thinner 1/2 oz
 PCS ISO 1/2 oz

STARTER KIT

Part num - 52-281
 6 Natural Shades in 1/2 oz bottle: PCS 2, PCS 3, PCS 5, PCS 7, PCS 9, PCS MILK
 PCS Thinner 1/2 oz
 PCS ISO 1/2 oz



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top Die Spacer



In order to control the relief and consistency of dies and castings, a uniform spacing before wax-up of a mold is essential. Our Die Spacer is a precise paint-on liquid that maintains a 12 micron homogeneous thickness, flows on smoothly and bonds to the gypsum or stone mold without running or clumping.

- | Precise control and viscosity to allow uniform coatings.
- | Fast drying, allowing quick application of multiple coats.
- | Contrasting color for definition and inspection.
- | Outstanding resistance to pressurized steam, chipping or process damage.

Part num	Color	Amount
26-9000	Blue	1 oz
26-900	Blue	1/2 oz
26-901	Red	1/2 oz
26-902	Gold	1/2 oz
26-903	Silver	1/2 oz
Thinner		
26-910	Transparent	1 oz

Top Die Lube

Top Die Lube is a fast-drying, non-oily, residue-free release film that is easy to apply to a die, providing a wet surface liner for wax castings. The thin finish allows for greater accuracy of the casting, no pattern distortion and clean separation, even after standing overnight. Excellent characteristics and performance off of stone, plated or epoxy dies.



- | Can be used with dripping wax or regular wax.
- | Easy removal, even after time, with no residue attached to casting or die.
- | Die spacer will not be removed with casting.
- | Contains alcohol, green color.

Part num	Color	Amount
26-800	Green	1 oz

Top Die Lube Pen Applicator

Pen Applicator -For the separation of plaster from wax

- | Reliable separating efficiency.
- | especially for the wax dipping technique.
- | Easy and quick removal of the wax-up from the model without distortions.
- | Economical Easy application due to pen shape.
- | Enables a uniform and thin application.
- | Adjusted to our Top waxes and Spacers.
- | Can be used with dripping wax or regular wax.

Part num	Color	Amount
26-801		1/2 oz
Refill		
26-802	Transparent	100ml



Top Rubbsep

Our Rubbsep is a high quality, easy to apply rubber film that forms a peel-able spacer between the die and the attachment during the molding process.



- | Provides accurate spacing for precision attachments.
- | Rubber coating prevents free monomer of resin from harming teeth.
- | Provides thermal protection for teeth.
- | Pick-up impression pour release agent for crowns.

Part num	Color	Amount
26-960	Pink	1/2 oz
26-961	White	1/2 oz

Top Blockout Wax



- | Fast and secure application to undercuts on dies for molding of comfortable attachments.
- | Adapts to die contours while soothing.
- | Does not burn without residue.
- | Will not coalesce with dipping wax, cervical wax or other modeling waxes due to very high melting point.
- | Does not discolor plaster model after boil-out.
- | For blocking out bite guard splints and dies.
- | 105 °C melting point.

Part num	Color	Amount
26-500	White Transparent	42/45g
26-501	Pink Transparent	42/45g







Ceramic Products

Top Ceraliq | Top Cera Glazer | Top Marker | Top Cera Sep
Top Die Hardener | Top Die Sealer | Top PC Spacer | Top PC Modeling Wax
Top Ceraflx Gel | Top Cera Contour | Ceramic Product Kits



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top Ceraliq



Our Top Ceraliq is a premier condensing liquid that assists any brand of porcelain powder to form a firmer, easier to control, hard finish. This product enhances the hardening properties, creating an ionic charge to automatically condense individual porcelain particles, helping them become denser and a tighter unit, enhancing the color, longevity and predictability of the application while eliminating air bubbles and particle clumping.

- | Longer moisture phase for predictable, reliable results.
- | No fear of dry-out, no need for frequent re-wetting and dilution.
- | Fast build-up, minimal shrinkage and more efficient condensing.
- | Shorter working time and a stronger, more natural looking result.

Part num	Color	Amount
26-510-16	Transparent	16 oz
26-510-32	Transparent	32 oz

Top Cera Glazer



Our Top Cera Glaze is a specially formulated liquid glaze medium that, when utilized during glaze firing results in a natural looking cosmetic restoration.

- | Universal Top Cera Glaze can be easily mixed with Hi Cera Solution for exact match.
- | Mix Top Cera Glaze liquid with glaze or stain powders for smooth consistency.

Part num	Color	Amount
26-520	Transparent	15 ml
26-521	Transparent	30 ml

Top Marker



A fast-drying low-viscosity marking pigment for detecting high spots and tight spots when seating castings, Top Marker provides high-contrast, accurate smears transferred at the slightest pressure on the opposing side to identify 'problem areas' on interproximal and occlusal areas.

- | Applies with brush.
- | Will not stain or discolor, does not contaminate metal.
- | Easy removal with water and soft brush.

Part num	Color	Amount
26-940	Blue	1/2 oz
26-941	Yellow	1/2 oz
Thinner		
26-945/1	Transparent	1/2 oz



Top Cerafix Gel

High quality holding gel for porcelain veneers & full-ceramic restorations on the die
 Top Cerafix Gel is a high quality temporary adhesive gel for holding crowns, Zirconia copings, veneers, wax and small metal castings in place during fabrication - procedures the build-up and contour processes.
 Carefully formulated with excellent hold veneers power, this gel is the perfect 'third hand' helper for temporary fitting.

- | Easy to Use
- | Top quality gel
- | Is especially useful in holding Pressed Veneers
- | Excellent holding and easy separation
- | Easily remove restoration from die
- | Burns without residue
- | Water-soluble
- | Transparent
- | Saves Time and Money

Part num	Color	Amount
26-515	Transparent	4 ml
26-516	Transparent	10 ml



Top Cera Contour



Top Cera Contour is a surface structures indicator. Ready for use by brush. Easily applied a thin layer. Visualization of imperfections, flaws, defects, irregularities allowing for the perfect match. Now you can achieve amazing results restorations in a short time and easily.

- | Defects and imperfection becomes visible
- | Simply and neatly smears
- | Dries within seconds
- | Easy clean with steam
- | Providing you the ability to create perfect restorations

Part num	Color	Amount
26-526	Gold	1/2 oz
26-527	Silver	1/2 oz



Top Die Hardener

Our Die Hardener is an absorbing resin that seals, hardens and strengthens stone models and dies. As it is absorbed into the stone surface, it permeates into all of the fissures and microscopic pores, bonding the entire stone by fusing the particles together. It resists chipping by leaving a zero micron infusion to the surface layer bonded to the resin permeating the stone, assuring an accurate and durable mold.



- | Appropriate for all stones and plasters, including gypsum.
- | Waterproofs and seals dental models and dies.
- | Resists chipping and cracking.
- | Reinforces margins to resist chipping.
- | Protects delicate detail and small features.

Part num	Color	Amount
26-930/4	Transparent	4 oz
26-930/8	Transparent	8 oz
26-930/16	Transparent	16 oz
Thinner		
26-911/4	Transparent	4 oz
26-911/8	Transparent	8 oz

Special thinner for setting preferred consistency Especially adapted to the texture of the TOP DIE HARDENER & SEALER

Top Die Hardener Pen Applicator

- | Reliable separating efficiency.
- | Comfortable and Easy application due to pen shape.
- | Economical and on-the-spot application;
- | Fast drying in only seconds. Hardens, seals, smooths
- | Fills open micro pores on the plaster surface
- | Abrasion-resistant surface. Dirt&water-repellent
- | Resists chipping and cracking.

Part num	Color	Amount
26-931	Transparent	1/2 oz
Refill		
26-932/4	Transparent	4 oz
26-932/8	Transparent	8 oz



Top Die Sealer

Our Top Die Sealer is brush applied to easily seal and water-proof stone die and model surfaces. It protects margins and reduces abrasion during operations, protecting the die from chipping, breaking or sticking. It also can provide an excellent base layer for our Top Cera Sep separating substance between porcelain and cast. It separates porcelain from model material on gingival surface of pontics, in both marginal porcelain regions as well as adjacent contact zones.



- | Protects preparation borders of crowns and nlays.
- | Prepares porcelain shoulder layers.
- | Adds negligible film thickness.
- | Thin consistency allows penetration into stone for firm seating.
- | Transparent surface becomes hard and water-proof.
- | Requires only one application, although additional thin layers can be applied after first coat is dry.

Part num	Color	Amount
26-920/4	Transparent	4 oz
26-920/8	Transparent	8 oz
Thinner		
26-911/4	Transparent	4 oz
26-911/8	Transparent	8 oz



Top Cera Sep

Our Top Cera Sep is perfect for sealing a die or model against the absorption of moisture when molding porcelain. Enhances release properties and prevents sticking of porcelain when separating from the die, prevents delicate, intricate parts from remaining with the die. Top Cera Sep does not cause discoloration of porcelain product, and burns out clean.



- | Separates porcelain from model material on gingival surface of pontics, in marginal porcelain regions and in adjacent contact zones.
- | Thin layers prevent end product faults.

Part num	Color	Amount
26-530	Transparent	1oz

Top Cera Sep Pen Applicator

Pen Applicator- For the separation of plaster from porcelain and wax
Recommended use on the sealed or hardened surface area

- | Reliable separating efficiency.
- | Easy application due to pen shape.
- | Easy and quick removal of the porcelain set-up
- | Prevents the moist porcelain to stick to plaster.
- | Enables a uniform and thin application.
- | Free of solvents.

Part num	Color	Amount
26-531	Transparent	1/2 oz

Reffil

26-532/4	Transparent	4 oz
26-532/8	Transparent	8 oz



Top Cera Sep Spray



Porcelain Release/Separating Agent.
Prevents porcelain from sticking to the stone model at the build-up.

- | Enhances release and prevents sticking of porcelain.
- | For release of onlays, occlusal & proximal contacts & resealing margins.
- | Conveniently packaged product - Pump Spray.
- | Easy to apply.
- | Does not cause discoloration.
- | Burns out clean.

Part num	Color	Amount
26-535/4	Transparent	4 oz
26-535/8	Transparent	8 oz



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top PC Spacer



Natural Tooth Colored- Stump Shade - Die Spacers System
FOR PRESSABLE CERAMIC TECHNIQUE

For the creation of highly aesthetic all-ceramic restorations, including inlays and onlays, veneers and partial crowns.

[More info in pages 9-10](#)

Part no	52-2801	52-2802	52-2803	52-2804	52-2805	52-2806	52-2807	52-2808	52-2809	52-2810
Color	PCS 1	PCS 2	PCS 3	PCS 4	PCS 5	PCS 6	PCS 7	PCS 8	PCS 9	PCS Milk
Amount	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz

Top PC Modeling Wax

Specifically formulated for press-able ceramic systems, our Top PC Modeling Wax is filtered several times to completely eliminate any opaque particles which might negatively react with the investment material or casted object. Like our Top Modeling Wax, this PC version retains excellent sculpting, caving and scraping characteristics due to our exclusive high density molecular structure.



- | No color contamination.
 - | Ash free for clean burnout.
 - | Compatible with all press-able systems.
 - | Low shrinkage.
 - | Fast application and quick cooling.
 - | High surface tension for precise carving and smooth finish.
 - | Compact cylinder contains 42-45g.
- Ash Free

Part no	26-400	26-401	26-402	26-403	26-404	26-405	26-406
Color	Grey	Blue	Green	Violet	Yellow	Bordeaux	Beige
Amount	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g



Ceramic Product Kits



Kit 1 Includes:

Part num - 168003

1. Top Cera Sep 1oz/Pen 0.5oz
2. Top Cera Glazer - 15ml
3. Top Ceraliq - 16oz

Kit 2 Includes:

Part num - 168004

1. Top Cera Sep Spray - 4oz
2. Top Cera Glazer - 15ml
3. Top Ceraliq - 16oz

Kit 3 Includes:

Part num - 168008

1. Top Marker Blue/Yellow - 0.5oz
2. Top Marker Thinner - 0.5oz
3. Top Die Sealer - 1oz

Top Ceraliq - universal modeling liquid to enhance superior porcelain results

Top Cera Glazer - improved cosmetic restoration at glance firing stage

Top Marker - high-contrast smears transferred at the slightest pressure

Top Cera Sep - better separation effect! prevents moisture absorption when build up porcelain



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com



DST – FLEXI DENTURE MATERIALS

DSTVITOlux-f | DSTVITOflexi | DSTVITOacry | DSTVITOcetal
DSTVITOkeep | DSTVITOpplast | DSTVITOCclasp | DST CLEAN & RENEW
DST LUSTERLIQ | DST ISOLIQ F | DST CEM F | DST FUSE FLEX | DST Lubriliq F



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DSTVITOflexi



Part no	52100	52101	52102	52110	52111	52112	52120	52121	52122
Color	Pink	Pink	Pink	Light Pink	Light Pink	Light Pink	Dark Pink	Dark Pink	Dark Pink
Amount	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL

Created from a unique thermoplastic, our **DSTVITOflexi** denture base material has exceptional aesthetic and physical properties, and is a rigid, yet flexible polymer ideal for creating compliant and comfortable partial dentures. It is moisture resistant, stain resistant and easy to adjust, finish and polish.

For the patient, selecting an appliance made from **DSTVITOflexi** results in a faster, more economical solution. Their partial denture will fit better, be more comfortable, unbreakable and easier to maintain.

For the dental technician, our polymer is easier to use, sets faster and the flexible properties mean accommodating a perfect fit with fewer complications. For dentists, prescribing **DSTVITOflexi** means the end result will be a stronger, flexible appliance which will be longer lasting and provide a superior function.

An exceptional alternative to metal framework fixtures. | Non-allergenic. | Non-brittle, virtually unbreakable. | High color stability in pink, light pink and dark pink. | Available in 25mm cartridges. | Cartridge size: medium, large, extra large. | 5 units per pack.



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DSTVITOacry



Part no	52200	52201	52202	52210	52211	52212	52220	52221	52222
Color	Pink	Pink	Pink	Light Pink	Light Pink	Light Pink	Dark Pink	Dark Pink	Dark Pink
Amount	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL

Our **DSTVITOacry** thermoplastic polymer is an advanced biocompatible material for injection-molding total and partial dentures. Featuring exceptional properties such as high strength, light weight, elasticity and natural translucence, it is more resistant against breakage than conventional resins. It does not cause allergenic reactions due to its monomer-free composition.

For the patient, selecting an appliance made from **DSTVITOacry** results in a faster, especially aesthetic and more economical solution. Their denture will fit better, be more comfortable and easier to maintain.

For the dental technician, our polymer is easier to use, sets faster and the aesthetic properties, which means accommodating a perfect fit with fewer complications.

For dentists, prescribing **DSTVITOacry** means the end result will be a stronger, flexible appliance which will be longer lasting and provide a superior function.

- Eliminates sensitivity to conventional acrylic dentures containing monomers. | Aesthetic, translucent and comfortable.
- Hypoallergenic and biocompatible. | Resistant to high impact. | High color stability in pink, light pink and dark pink.
- Available in 25 mm cartridges. | Cartridge size: medium, large, extra large. | 5 units per pack.



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DSTVITOCetal



Part no	52300	52301	52302	52310	52311	52312	52320	52321	52322
Color	Pink	Pink	Pink	A2	A2	A2	A3	A3	A3
Amount	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL

DSTVITOCetal is a biocompatible, non-allergenic and pure, crystalline acetyl copolymer resin which is used for partial and total denture frameworks without the need for cast metal backbones. With very low moisture absorption properties, an appliance cast from this material is very low weight, maintains high tensile strength and excellent flexural strength, hardness and fatigue resistance.

Can also be utilized for the creation of clasps, provisional crowns and bridges, and unilaterals.

For the patient, selecting an appliance made from **DSTVITOCetal** results in a faster, more economical solution. Their partial denture will fit better, be more comfortable, unbreakable and easier to maintain.

For the dental technician, our polymer is easier to use, sets faster and the flexible properties mean accommodating a perfect fit with fewer complications.

For dentists, prescribing **DSTVITOCetal** means the end result will be a stronger yet lighter and more flexible appliance which will be longer lasting and provide a superior function.

- Aesthetic and comfortable. | Superb dimensional stability and elastic memory. | Low moisture absorption. | 5 units per pack.
- High color stability in pink, A2, A3 | Available in 25 mm cartridges. | Cartridge size: medium, large, extra large.



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DSTVITOkeep



Part no	52400	52401	52440	52441
Color	Tan	Tan	Pink	Pink
Amount	50gr	100gr	50gr	100gr

DSTVITOkeep is a semi-crystalline thermoplastic with excellent mechanical and chemical-resistance properties. As a high performance polymer is the best alternative for non-metallic restorations to fabricate full and partial removable dentures and implant overdentures frameworks, crown and bridge construction.

- Novel biomaterial as an alternative to zirconia
- Individual adjustment and prevents chipping
- Optimal work flow
- Allows a high standard of precision and finishing
- Perfect aesthetic results
- Magnificent feeling
- High resistance to abrasion
- Ideal mechanical properties
- Particularly attractive
- Biocompatibility, Greater Versatility



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DSTVITOp^last

Softens in hot water

NEW



Part no	52700	52701	52702	52710	52711	52712
Color	Blue	Blue	Blue	Yellow	Yellow	Yellow
Amount	300 gr	500 gr	1 kg	300 gr	500 gr	1 kg
Part no	52720	52721	52722	52730	52731	52732
Color	Pink	Pink	Pink	Green	Green	Green
Amount	300 gr	500 gr	1 kg	300 gr	500 gr	1 kg

DSTVITOp^last is an innovative alternative for conventional custom tray material. Designed for a variety bite registration -fast and easy to use. Friendly and convenient for the patient and dentist - comfort custom trays, bite registration, fore and side bite detection in 1.2.3. Features good working time along with excellent malleability and adaptation to model. Easy trimming characteristics reduce effort.

- Affordable, cost effective
- Biocompatible and non-allergenic
- Newest alternative for conventional tray material
- Softens in hot water
- Rigid yet unbreakable
- Comfortably handled and shape
- Can be simply cut and fit in softened state
- Extraordinary accuracy
- No odor, very durable



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DSTVITOlux-f



Part no	52010	52011	52012	52000	52001	52002	52020	52021	52022
Color	Light Pink	Light Pink	Light Pink	Pink	Pink	Pink	Dark Pink	Dark Pink	Dark Pink
Amount	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL	25mm M	25mm L	25mm XL

This innovative flexible dental resin allows users to easily create comfortable and compliant removable dentures. Created from a unique thermoplastic material, DSTVITOlux-f has exceptional aesthetic and physical properties that enable highly comfortable and natural looking partials. This rigid, yet flexible polymer is an ideal solution.

DSTVITOlux-f's biocompatible, high-performance polymer offers all the advantages of acrylic, enabling the creation of aesthetic, natural-looking comfortable prosthetic applications, including full and partial dentures, splints, clasps and more. Resistant to both moisture and stains, DSTVITOlux-f is easy to adjust, finish and polish. It sets faster than other materials, and because of its flexible properties, patients enjoy an improved fit.

- Revolutionary, Rigid yet flexible denture base material
- Easily repaired and relined by conventional acrylic
- luxury comfort, durability and moderate flexibility
- Strong yet flexible polymer is easy to work with
- Repels plaque and odor
- Exceptional aesthetic and physical properties
- Enables a more perfect fit



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

DST CLEAN & RENEW

The DST Clean/Renew concentrated denture cleaner powder, is an excellent solution and recommended for cleaning thermoplastic and acrylic dentures.



- | Do not apply while denture is in mouth, use as instructed with dentures removed.
- | Can be mixed and used for up to 7 days before being replaced with a new mix.
- | One spoonful of powder is mixed with 7 oz. of warm water.

Part num		Amount
52-690	powder	4 oz

DST LUSTERLIQ

DST LUSTERLIQ can be used to the finish gloss on any acrylic and thermoplastic denture appliances.



Part num	Color	Amount
52-620/4	Transparent	4 oz
52-620/8	Transparent	8 oz

- | Easy to apply
- | resulting in the quickly, high-gloss luster.
- | Provides a brilliant surface finish
- | Leaves a pleasant fresh taste

DST LUBRILIQ F

DST Lubriliq it is perfect for lubricating aluminum cartridges prior to insertion into a press oven, or for greasing the flasks.



Part num	Color	Amount
52-610	Transparent	8 oz

- | Multi-purpose release lubricant
- | easy and safe release
- | Suitable for all injection molding machines
- | Easy to application assures even coverage



DST ISOLIQ F



Isolating and separating the plaster mold from the material is important in protecting both the die from damage and the appliance from debris and powder castoff. Appliance separates in a shiny condition, free from plaster debris and requiring minimum finishing.

The use of DST Isoliq as an isolation compound that offers several advantages, including water repellent characteristics, a protective layer which isolates the gypsum from thermoplastics, and high heat resistance, which makes it ideal for molding injection processes which rely on heat and pressure.

- | Seals and repels water due to its thin isolation layer.
- | Resistant to high temperature.
- | Fast drying liquid applies with brush

Part num	Color	Amount
52-600	Transparent	8 oz
52-601	Transparent	16 oz

DST CEM F



Our DST Cem F is a very strong transparent glue that bonds artificial acrylic teeth into the stone or the flask securely and firmly, after boil out - prior to injection. DST Cem F is mild enough to keep from bonding the teeth to the stone post-injection.

- | Easy use with applicator of thin liquid glue.
- | Fast drying without residue.
- | Will not stain or damage surface of teeth.

Part num	Amount
52-650	20 ml

DST FUSE FLEX



This is a liquid fusing bond that is formulated for repairs, linings or additions to polyimide resin appliances and partials. DST Fuse Flex begins to deteriorate existing bonds, and creates a new bond when injected into the area. The new bond that is created is stronger and virtually undetectable.

- | Makes repairs and additions to existing partials immediately prior to re-injection.
- | Use for minor repairs and additions to existing DST Unbreakable Partial
- | Use only when flasking and injection are absolutely necessary.
- | Not to be used in the mouth. Must be used in dental laboratory

Part num	Amount
52-640	1 oz





Denture Products

Top Natur Glaze | VITO-CRYL Seal - IQ | Top Airstain Glaze
Top Air Opaque | Top Paqubond | Top Pumiliq | Top Gyps-Rid | Top Isoliq A
Top Blockout Wax | Top Bite Rims | Top Sticky Wax | Top Spec Sticky Wax



Dental Materials Production
www.dst-dental.com | Email: info@dst-dental.com



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top Natur Glaze

Clear acrylic liquid for surface sealing of a denture and any acrylic appliance .

Brilliant natural finish, without mechanical polishing.

Top Natur Glaze produces a natural wet-looking protective glaze over permanent and provisional acrylic allowing the base color to show through. Self-curing, increases color stability. Fills micro voids to reduce bacteria and odor problems.

Bonds to all dental resins, including methyl and ethyl methacrylates, polycarbonates, bis-acryl resins, resin crowns and laminates, vinyl methacrylate copolymers, acrylic denture teeth denture base, and composite resins.



- | Brilliant natural finish, without mechanical polishing.
- | Fully transparent, no additional polishing required.
- | Creates a very smooth surface.
- | High gloss, colourless, strong bonding
- | Also protects against plaque deposits and discolouration.
- | Applicable for composite restoration.
- | Rapid air dry within 20 seconds.

Part num	Color	Amount
26-272	Transparent	1/2 oz
Thinner		
26-275	Transparent	1/2 oz

VITO-CRYL Seal - IQ®

Clear acrylic spray to seal & shine the tissue surface a denture and any acrylic appliance

Clear acrylic spray to seal&shine a fills open pores the tissue side of a denture and any acrylic appliance, to seal the tissue surfaces from both bacteria and plaque and Brilliant natural finish, without rough spots.



- | EASY TO USE!
- | Impressive results in less time
- | Micro thin layer
- | Creates an excellent Seal&Shine smooth surface
- | Improve patient comfort,
- | Dries quickly - 60 to 90 seconds
- | Chemically bonds to all-types of acrylic resin

Part num	Color	Amount
18-101	Transparent	1 oz
18-102	Transparent	2 oz



Top Air Opaque



Our Top Air Opaque is a colored opaque liquid which performs a secure bond between metal surfaces/framework and acrylic or other resins. It is brushed on to achieve desired thickness and air dries quickly. Its pink tint renders it almost invisible when the opaque is applied to the framework.

The bonding properties of Top Air Opaque can be enhanced by mixing in Top Paqubond, a powder that is mixed with the Top Air Opaque to create a harder film which will not dislodge during acrylic packing, and the resulting bonding film will then be resistant to softening by the monomer, eliminating dislodging during the process. This compound can be used under cured resins, acrylics, bis-acrylic and composite resins.

Part no	26-2901	26-2902	26-2903	26-2904	26-2905	26-2906
Color	Pink	Light Pink	Dark Pink	Snow	Cream - A+B	Beige - C+D
Amount	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz	1/2 oz

Top Air Opaque Thinner - Part no -26-295, Color - Transparent, Amount - 1/2 oz

Top Paqubond



Top Paqubond is a specially formulated opaquer basic agent - comprises powder&liquid system used to enhance the bond to TOP AIR OPAQUE, and outstanding adhesion values in the case of compressive, shearing and tensile loads for use in dental restorations where it is necessary to mask the appearance of metal for pleasing aesthetics.

	POWDER	LIQUID	SET
Part no	26-2922	26-2921	26-2923
Amount	7g	10ml	

Top Air Opaque - KITS



COMPLETE KIT

Part num - 26-290
6 colors in 1/2 oz Bottle - Pink, Light Pink , Dark Pink, Snow, Cream - A+B, Beige - C+D
Top Air Opaque Thinner - Transparent - 1/2 oz
Top Paqubond- Special to Top Air Opaque - POWDER 7g
Top Paqubond- Special to Top Air Opaque - LIQUID 10ml
Mixing Pad
Micro Applicator Brushes

STARTER KIT

Part num - 26-291
3 colors in 1/2 oz Bottle - Pink, Snow, Cream - A+B,
Top Air Opaque Thinner - Transparent - 1/2 oz
Top Paqubond- Special to Top Air Opaque - POWDER 7g
Top Paqubond- Special to Top Air Opaque - LIQUID 10ml
Mixing Pad
Micro Applicator Brushes



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top Airstain Glaze



Intended for esthetic color modification of any acrylic appliance.
Effects of natural realistic appearance.

Our Top Air Stain Glaze is an acrylic stain, creating a natural finish for quick characterizing, shading and staining of provisional and permanent acrylic and composite restorations - denture base, acrylic denture teeth and composite resins, for a more realistic appearance.

Top Air Stain Glaze are available in a 11 natural colors for perfect blending and matching;
 All 11 Art stains can be applied alone or blended together to achieve different effects, available for color characterization and the shade effects of natural teeth and the pigmentation of the gums.

Translucent and stable, these stains do not change tint or hue during application or over time. It cures quickly in less than 20 seconds, and the color can be built up gradually.

for a more realistic appearance | Excellent properties bonds and cure to all dental resins | applied by brush
Sealing properties fill porosit | Also available in a choice Kit

Available
 In choice 9 individual colors

Part no	26-2700	26-2701	26-2702	26-2703	26-2704	26-2705
Color	Violet	White	Brown	Gray	Blue	Pink
Amount	10ml	10ml	10ml	10ml	10ml	10ml
Part no	26-2706	26-2707	26-2708	26-2709	26-2710	
Color	Dark Pink	Yellow	Vanilla- A+B	Beige- C+D	Orange	
Amount	10ml	10ml	10ml	10ml	10ml	

Complementary products

	Top Airstain/Natur Glaze Thinner	Thinner Brush Cleaning
Part no	26-275	26-274
Color	Transparent	Transparent
Amount	1/2 oz	1/2 oz

Top Airstain Glaze - KITS



Professional Kit -9 Colors

Part num - 26- 270

9 colors in 10ml Bottle - Pink, Dark Pink, Yellow, Brown, Gray, Blue, White and two cervical colors: Vanilla- A+B, Beige- C+D

Top Natur Glaze - clear - 1/2 oz

Thinner - 1/2 oz

Brush Cleaning Thinner - 1/2 oz

1 Mixing Well

2 Brush

Tooth Aesthetic Kit - 4 Colors

Part num - 26- 273

4 colors in 10ml Bottle - Gray, Yellow, Vanilla- A+B, Beige- C+D

Top Natur Glaze - clear - 1/2 oz

Thinner - 1/2 oz

Brush Cleaning Thinner - 1/2 oz

1 Mixing Well

2 Brush

Denture Aesthetic Kit - 4 Colors

Part num - 26- 271

4 colors in 10ml Bottle - Pink, Dark Pink, Blue, White

Top Natur Glaze - clear - 1/2 oz

Thinner - 1/2 oz

Brush Cleaning Thinner - 1/2 oz

1 Mixing Well

2 Brush



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top Pumiliq



Our Top Pumiliq is a mixing liquid for polishing compounds, containing a disinfectant and used with any and all pumice powders. It maintains a creamy consistency for a full day, streamlining work processes. Top Pumiliq contains small amounts of stabilizing agents.

- | Skin protective.
- | Free from formaldehyde.
- | Stays in creamy working consistency all day.

Part num	Color	Amount
51-220-1000	Transparent	32 oz

Top Gyps-Rid



- | Highly effective gypsum remover in powder form.
- | Becomes fully efficient in a 35-45 °C ultrasonic bath.
- | Stirring the powder with water creates a solution to the viscosity required.
- | Dilution ratio is 100g powder to 400-500ml water.
- | Keeps the utensils shiny and clean.
- | Active material concentration can be customized by powder/water proportions.

Part num	Color	Amount
26-500	Powder	500 g

Top Blockout Wax

A very hard, high temperature wax, our Top Blockout Wax is designed to fill undercuts and process induced cavities, cover sharp edges and smooth out rough surfaces of dies prior to molding attachments. Excellent carvability and very good scraping properties.



- | Fast and secure application to undercuts on dies for molding of comfortable attachments.
- | Adapts to die contours while soothing.
- | Does not burn without residue.
- | Will not coalesce with dipping wax, cervical wax or other modeling waxes due to very high melting point.
- | Does not discolor plaster model after boil-out.
- | For blocking out bite guard splints and dies.
- | 105 °C melting point.

Part num	Color	Amount
26-511	Light Pink- opaque	42/45g
26-501	Pink Transparent	42/45g

Top Bite Rims

Anatomically formed, flexible and comfortable, our Top Bite Rims Wax is able to achieve exact bite registration in one application. These pre-formed rims include a basal channel for good fit and registration, and come in a hardness that accepts excellent tooth detail, yet is comfortable to bite.



Part num	Color	Amount
26-130	Red-flavored preformed	50 pcs

- | Flavored and comfortable for patient.
- | Standard jaw sizes.
- | Contains 50 pieces per box.

Top Isiliq A



PLASTER & ACRYLIC SEPARATOR

Tin foil substitute for separation of acrylic appliances from stone models

TOP ISOLIQ A is an alginate-based separating solution formulated for use between gypsum stone & plaster and acrylic resin. For the insulation of plaster surfaces against denture resins (heat and cold curing) and composites.

- | Prevents cured acrylic denture from sticking to stone plaster cast.
- | Thin separating film.
- | Short action time and no formation of layers
- | Achieving smooth surfaces
- | Optimum separation - High quality
- | Works with cold or hot models

Part num	Color	Amount
26-750/16	Pink	16oz
26-750/32	Pink	32oz

Top Sticky Wax



Adhesive power is one of the remarkable properties of our Top Sticky Wax -making it the best choice for break repairs, fixing teeth, front walls and solder molds. Quick setting and high strength prevent distortion and slippage during work operations. It is essentially a 'wax glue', very tacky with low flexibility after setting, making it a top choice for positioning, holding in place and making sprue connections.

- | 'Wax glue' properties for positioning.
- | Stays viscous to make attachments.
- | High stability after hardening maintains dimensional integrity.
- | Can be cleaned up or removed with steam jet.
- | Leaves no residue or contamination.

Part num	Color	Amount
26-250	Neon-Pink	42-45g

Top Spec Sticky Wax



Top Spec Sticky Wax with extremely strong adhesive properties, is formulated to meet the needs of professionals. It quickly and accurately fixates acrylic teeth within a silicone mould. This product was explicitly developed by our specialists and allows the technician to fixate acrylic teeth without the use of glues, assuring perfect results.

- | Fixates acrylic teeth - no adhesive application.
- | Prevents damage to the anatomical shape of acrylic teeth.
- | Quick, easy and safe to apply.
- | Extreme hardness, distortion is impossible.
- | Easily removed with steamers.
- | Leaves no residue.

Part num	Color	Amount
26-300	Red	42-45g





Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com



Cad/Cam

Top Scan Wax | Bari x-ray



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Top Scan Wax

CAD/CAM opaque

Our Top Scan Wax for CAD/CAM is a specially formulated laser-opaque and light-opaque wax for sculpturing according to the CAD/CAM procedure. Guarantees superb scanning results which enhances exact fitment. It builds up fast, sets rapidly and is quickly available for carving and scraping operations.



- | Highly opaque for prime scan operation.
- | No powder spray necessary.
- | High stability after hardening maintains dimensional integrity.
- | Superb modeling qualities.
- | Not to be used in casting.

Part num	Color	Amount
26-550	Beige- opaque	42-45 g
26-551	purple - opaque	42-45 g
26-552	Turquoise- opaque	42-45 g

Bari x-ray

To mix with acrylics for X-ray feature, for optimum diagnostic imaging.

Bari x-ray is excellent for adding and mixing with Acrylic Resin, provides a high contrast identifier in X-Ray scanning pictures with consistent texture, for the fabrication of methacrylate - based acrylic appliances and used wherever a radiopaque acrylic dental device is useful in establishing an anatomical location or site , such as a provisional implant stent.



- | Ideal for production of radiopaque acrylic dental device
- | diagnostics, and planning of treatments.
- | Identifies primary implant location.
- | borders are sharp for precise implant planning

Part num	Color	Amount
52-700	Powder	4 oz



Wax Products

Dipping Waxes | Crown and Bridge Waxes | Special Waxes
Bite Registration Waxes | Sticky Waxes



Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com

Dipping Waxes

Top Dip Wax is used in electrical dipping pots. Perfect measures allow precise and consistent results. Once liquefied, provides excellent elasticity, adaption and detail. Retains strong edges with a memory enhanced elastic feel, clean burn out and minimum shrinkage for consistent dimensional integrity.

- | Dries quickly.
- | Eliminates resealing.
- | Maximum hygroscopic expansion for casting fit requiring no additional adjustment.
- | Contains natural resins for stability and elasticity.
- | Excellent scraping and cutting properties allow best fittings and easy crown adaption.
- | Resists warping, deforming or breaking of copings.
- | Can be held in molten state for long periods of time.
- | Easy clean up.

Top Dip Wax S



Part num	Color	Amount
26-100D	Red	100g

- | Homogeneous wax form model results.
- | Leaves no wax residue on model.
- | Retains strong edges.

Top Dip Wax M



Part num	Color	Amount
26-200D	Orange	100g
26-210D	Pink	100g

- | Perfect elasticity properties.
- | Recommended hardness for all applications.

Top Dip Wax H



Part num	Color	Amount
26-300D	Green	100g

- | Great for making secondary forms with conus or telescope technique.
- | Precise and consistent results.

Crown and Bridge Waxes

Top Modeling Wax

Our Top Modeling Wax features optimum properties for carving and scraping, due to non-elasticity, hardness without brittleness, and high surface tension. This wax has quick application and cooling properties, speeding up the work flow when building bridges, crowns and inlays.



- | Optimal raw materials allow additive wax-up techniques.
- | High density molecular structure is easily carvable and leaves smooth surfaces.
- | Strong enough for elaborating detailed cusp tips.
- | Low shrinkage, low expansion, opaque.
- | Clean burnout, ash free, leaves no traces.

Part no	26-200	26-201	26-202	26-203	26-204	26-205	26-206
Color	Gray- opaque	Blue- opaque	Green- opaque	Violet- opaque	Brown- opaque	Red- opaque	Beige - opaque
Amount	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g

Top PC Modeling Wax

Specifically formulated for press-able ceramic systems, our Top PC Modeling Wax is filtered several times to completely eliminate any opaque particles which might negatively react with the investment material or casted object. Like our Top Modeling Wax, this PC version retains excellent sculpting, caving and scraping characteristics due to our exclusive high density molecular structure.



- | No color contamination.
- | Ash free for clean burnout.
- | Compatible with all press-able systems.
- | Low shrinkage.
- | Fast application and quick cooling.
- | High surface tension for precise carving and smooth finish.
- | Ash Free

Part no	26-400	26-401	26-402	26-403	26-404	26-405	26-406
Color	Gray	Blue	Green	Violet	Yellow	Bordeaux	Beige
Amount	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g

Top Sculpturing Wax

Offering a very firm consistency that allows excellent handling and carving suitability, our Top Sculpturing Wax will not become brittle, crack or fracture during process activities or curing. This is a harder wax that is easy to scrape, carve and create detailed occlusal anatomy. When heated and in pouring or dipping form, the viscosity is easily controlled, covering large areas expediently and efficiently.



- | Low contraction factor allows precision fit.
- | Low shrinkage for unsurpassed fabrication accuracy.
- | Clean burnout due to resin and ash free manufacture.
- | Stable inner elasticity allows for true reproduction of occlusal stops and contact points.
- | Ash Free

Part no	26-100	26-101	26-102	26-103	26-104	26-105	26-106
Color	Gray	Blue	Neon Green	Neon Yellow	Orange	White	Beige
Amount	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g	42-45 g

Special Waxes

Top Diagnostic Wax

Specially formulated for diagnostic use, our Top Diagnostic Wax has an elevated melting point and a natural finish for lifelike results. Modeling properties match the best properties of our leading sculptural waxes while providing extra elasticity for accurate representation of detail, surface and finish.



Part num	Color	Amount
26-450	Dentine-opaque	42-45 g
26-451	Write-transparent	42-45 g
26-452	Write-opaque	42-45 g

- | Extremely stable, low in contraction.
- | Ideal for elaborating wax-ups.
- | Easy to scrape and carve.
- | Color corresponds to tooth enamel.
- | Not to be used in casting.

Top Cervical Wax

Ensuring a perfect marginal seal for inlays, onlays and crowns requires our completely shrink-free Top Cervical Wax. Adheres perfectly and comfortably to the coping wax, its smooth flow provides precise margins and secure fitment. Carving and scraping touch up can be carried out without tearing, allowing construction of very thin margins.



Part num	Color	Amount
26-350	Red	42-45 g

- | Non-rebounding and tension-free.
- | No pulling or flaring.
- | Compatible with Top Dip Wax and Top Sculpturing Wax.
- | Easy flowing, non-smear qualities allow precision application and a sharp and accurate margin.
- | Ash free.



Special Waxes

Top Milling Wax

Our Top Milling Wax is an extra hard wax that is configured with sculptural handling characteristics as well as the best properties of a high quality milling wax. Designed for customizing implant components and abutments as well as every aspect of ceramic crown and bridge waxing, it can also be used for pressed ceramics due to clean burn out.



- | Superb, stable modeling properties.
- | Ash and Resin free for clean burnout properties.
- | Stable and tough to accept variable speed milling tools.
- | No wax smudging, smearing or sticking to tools during burring.
- | Burs can create detail; the result will be smooth and clean.

Part num	Color	Amount
26-150	Gray - opaque	42-45 g
26-151	Blue - Transparent	42-45 g
26-152	Gray- Transparent	42-45 g

Top Scan Wax

Our Top Scan Wax for CAD/CAM is a specially formulated laser-opaque and light-opaque wax for sculpturing according to the CAD/CAM procedure. Guarantees superb scanning results which enhances exact fitment. It builds up fast, sets rapidly and is quickly available for carving and scraping operations.



- | Highly opaque for prime scan operation.
- | No powder spray necessary.
- | High stability after hardening maintains dimensional integrity.
- | Superb modeling qualities.
- | Not to be used in casting.

Part num	Color	Amount
26-550	Beige- opaque	42-45 g
26-551	purple - opaque	42-45 g
26-552	Turquoise- opaque	42-45 g

Top Blockout Wax

A very hard, high temperature wax, our Top Blockout Wax is designed to fill undercuts and process induced cavities, cover sharp edges and smooth out rough surfaces of dies prior to molding attachments. Excellent carvability and very good scraping properties.



- | Fast and secure application to undercuts on dies for molding of comfortable attachments.
- | Adapts to die contours while soothing.
- | Does not burn without residue.
- | Will not coalesce with dipping wax, cervical wax or other modeling waxes due to very high melting point.
- | Does not discolor plaster model after boil-out.
- | For blocking out bite guard splints and dies.
- | 105 °C melting point.

Part num	Color	Amount
26-511	Light Pink- opaque	42-45 g
26-501	Pink Transparent	42-45 g
26-500	white opaque	42-45 g
26-510	white Transparent	42-45 g

Bite Registration Waxes

Top Bite Rims



Part num	Color	Amount
26-500	Red-flavored preformed	50 pcs

- | Flavored and comfortable for patient.
- | Standard jaw sizes.
- | Contains 50 pieces per box.

Sticky Waxes

Top Sticky Wax

Adhesive power is one of the remarkable properties of our Top Sticky Wax -making it the best choice for break repairs, fixing teeth, front walls and solder molds. Quick setting and high strength prevent distortion and slippage during work operations. It is essentially a 'wax glue', very tacky with low flexibility after setting, making it a top choice for positioning, holding in place and making sprue connections.



- | 'Wax glue' properties for positioning.
- | Stays viscous to make attachments.
- | High stability after hardening maintains dimensional integrity.
- | Can be cleaned up or removed with steam jet.
- | Leaves no residue or contamination.

Part num	Color	Amount
26-250	Neon-Pink	42-45g

Top Spec Sticky Wax

Top Spec Sticky Wax with extremely strong adhesive properties, is formulated to meet the needs of professionals. It quickly and accurately fixates acrylic teeth within a silicone mould. This product was explicitly developed by our specialists and allows the technician to fixate acrylic teeth without the use of glues, assuring perfect results.



- | Fixates acrylic teeth - no adhesive application.
- | Prevents damage to the anatomical shape of acrylic teeth.
- | Quick, easy and safe to apply.
- | Extreme hardness, distortion is impossible.
- | Easily removed with steamers.
- | Leaves no residue.

Part num	Color	Amount
26-300	Red	42-45g



CASTING/ SOLDERING

Top Wax Finish | Uni-Wax Adhesive | Top Vacudeb



Top Wax Finish

Our Wax Finish is a premium smoothing wax sealer, applied to the occlusal surface of a modellation; this product smooths out potential irregularities and levels small scratches and unevenness. Floods any fissure depths that can be inaccessible for a fissure frill.



- | Creates a very smooth, monogeneous surface.
- | Dries in 2-3 minutes, depending on thickness of layer.
- | Finer surface after modelling reduces elaboration time by up to 50%.
- | Occlusal and approximal contact points remain constant after polishing.
- | Reinforcement of contact points.

Part num	Color	Amount
26-000	Transparent	1/2 oz
Thinner		
26-005	Transparent	1/2 oz

Uni-Wax Adhesive

For use with preformed wax patterns, our Uni Wax Adhesive is a residue free bonding agent for use with preformed wax patterns used with unhardened investment duplicate models. Adhesive can be applied thinly to ensure retention of patterns on the duplicate model or the investment model without any gap or airspace, providing a secure adherence. Thin application allows part of the liquid to soak into the investment material to strengthen the bond.



- | Useful for wax and resin components.
- | Contains no heavy metals.
- | Leaves no film or residue.
- | Burns without residue.

Part num	Color	Amount
26-050	Transparent	1/2 oz

Top Vacudeb

A surface tension reducer, Top Vacudeb allows an attachment to achieve close adaption and maximum contact with all surfaces of a wax pattern and with all impressions prior to investing, without trapping air, forming air bubbles or voids on the surface of models and castings.



- | May be used under vacuum.
- | Fast drying increases productivity.
- | Eliminates air entrapments, bubbles and voids.
- | Facilitates fluid flow to produce sharper detail.
- | Supports silicone, wax, rubber and alginate patterns.

Part num	Color	Amount
26-220-250	Pink- Transparent	8 oz
26-220-500	Pink- Transparent	16 oz
26-220-1000	Pink- Transparent	32 oz





Dental Materials Production

www.dst-dental.com | Email: info@dst-dental.com



DS TECHNOLOGY LTD

29 Omarim st., Omer Industrial Park, 8496500, Israel
info@dsdentalisrael.com | www.dsdentalisrael.com
Tel. +972-77-2006636 Fax. +972-72-2829896